


## 200 MHz Variable Gain Photoreceiver




Features

- Adjustable transimpedance gain from  $10^2$  to  $10^8$  V/A
- Wide bandwidth up to 200 MHz
- Si-PIN photodiode covering the 400 to 1000 nm wavelength range
- Large optical detector size 1 × 1 mm
- High dynamic input range up to 10 mW optical power
- Very low noise, NEP down to 76 fW/√Hz
- Switchable low pass filters for minimizing wideband noise
- Free-space input 1.035"-40 threaded, easily convertible to fiber optic input (FC and FSMA) with optionally available screw-on adapters
- Full manual and remote control capability

Applications

- All-purpose low-noise photoreceiver (O/E converter) for the MHz range
- Time resolved optical pulse and power measurements
- Laser intensity noise measurements (RIN)
- Optical front-end for oscilloscopes, spectrum analyzers, A/D converters and RF lock-in amplifiers

Block Diagram



BS01-OE-300\_R2

## 200 MHz Variable Gain Photoreceiver

**Intended Use**

The OE-300-SI-10 is a high speed variable gain photoreceiver. It is designed for fast and precise conversion of small optical signals into equivalent output voltages. Operation is mostly self-explanatory. If in doubt, consult this document or contact support@femto.de.

For safe operation, please refer to the damage thresholds specified in the "Absolute Maximum Ratings", "Temperature Range" and "Power Supply" sections of this document.

The operating environment must be free of smoke, dust, grease, oil, condensing moisture, and other contaminants that could affect the operation or performance.

**Available Version**

OE-300-SI-10-FST



1.035"-40 threaded flange with internally threaded coupler ring (outer diameter 30 mm) for free space applications. Compatible with many optical standard accessories and for use with various types of fiber connector adapters.

Optionally available:

Fiber adapters PRA-FC, PRA-FCA and PRA-FSMA. With the relative large 1 x 1 mm photodiode installed in the OE-300-SI-10 input coupling is not critical. However, standard SM 9/125 fibers (PC or APC) with low numerical aperture (NA) are recommended for ensuring near 100% coupling efficiency.

**Related OE-300 Models**

See separate datasheets for following models on [www.femto.de](http://www.femto.de):

OE-300-SI-30-FST

Si-PIN, Ø 3 mm, 320 - 1000 nm  
1.035"-40 threaded flange

OE-300-IN-01-FC

InGaAs-PIN, Ø 80 µm, 900 - 1700 nm  
FC fiber receptacle only

OE-300-IN-03-FST

InGaAs-PIN, Ø 300 µm, 800 - 1700 nm  
1.035"-40 threaded flange

**Available Accessories**

PRA-FC  
PRA-FCA  
PRA-FSMA



Fiber-adapter with external 1.035"-40 thread

PRA-PAP




Alternative mounting option: post adapter plate, easy to mount on FEMTO photoreceiver series OE, FWPR, PWPR, HCA-S and LCA-S

PS-15-25-L




Power Supply  
input: 100 – 240 VAC  
output: ±15 VDC

## 200 MHz Variable Gain Photoreceiver

Available Accessories (continued)	LUCI-10	<div style="text-align: center;">  </div> <p>Compact digital I/O interface for USB remote control, supports opto-isolation of amplifier signal path from PC USB port, 16 digital outputs, 3 opto-isolated digital inputs, bus-powered operation</p>																																																																																																																																																	
Specifications	<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 30%;">Test conditions</td> <td colspan="6"><math>V_S = \pm 15\text{ V}</math>, <math>T_A = 25\text{ }^\circ\text{C}</math>, output load impedance <math>50\ \Omega</math>, warm-up 20 minutes (min. 10 minutes recommended)</td> </tr> <tr> <td rowspan="2">Gain</td> <td>Transimpedance gain</td> <td colspan="5"><math>1 \times 10^2 \dots 1 \times 10^8\ \text{V/A}</math> (output load <math>50\ \Omega</math>)</td> </tr> <tr> <td>Gain accuracy</td> <td colspan="5"><math>\pm 1\%</math> electrical, between settings</td> </tr> <tr> <td rowspan="2">Frequency Response</td> <td>Lower cut-off frequency</td> <td colspan="5">DC / 100 Hz, switchable</td> </tr> <tr> <td>Upper cut-off frequency (-3 dB)</td> <td colspan="5">up to 200 MHz (see table below), switchable to 1 MHz or 10 MHz</td> </tr> <tr> <td rowspan="2">Input</td> <td>Optical CW saturation power</td> <td colspan="5">see table below</td> </tr> <tr> <td>Noise equivalent power (NEP)</td> <td colspan="5">see table below</td> </tr> <tr> <td rowspan="20">Performance depending on Gain Setting</td> <td><u>Gain setting (low noise) (V/A)</u></td> <td><math>10^2</math></td> <td><math>10^3</math></td> <td><math>10^4</math></td> <td><math>10^5</math></td> <td><math>10^6</math></td> <td><math>10^7</math></td> </tr> <tr> <td>Upper cut-off frequency (-3 dB)</td> <td>200 MHz</td> <td>80 MHz</td> <td>14 MHz</td> <td>3.5 MHz</td> <td>1.8 MHz</td> <td>220 kHz</td> </tr> <tr> <td>Rise/fall time (10 % - 90 %)</td> <td>1.9 ns</td> <td>3.25 ns</td> <td>26 ns</td> <td>92 ns</td> <td>235 ns</td> <td>1.6 <math>\mu\text{s}</math></td> </tr> <tr> <td>NEP (<math>\sqrt{\text{V/Hz}}</math>, @850 nm)</td> <td>322 pW</td> <td>25 pW</td> <td>2.9 pW</td> <td>740 fW</td> <td>260 fW</td> <td>78 fW</td> </tr> <tr> <td>Measured at</td> <td>20 MHz</td> <td>8 MHz</td> <td>1.4 MHz</td> <td>350 kHz</td> <td>180 kHz</td> <td>22 kHz</td> </tr> <tr> <td>Integr. input noise (RMS)*</td> <td>7.5 <math>\mu\text{W}</math></td> <td>580 nW</td> <td>35 nW</td> <td>4.9 nW</td> <td>1.3 nW</td> <td>100 pW</td> </tr> <tr> <td>CW saturation power (@ 850 nm)</td> <td>10 mW</td> <td>1.7 mW</td> <td>170 <math>\mu\text{W}</math></td> <td>17 <math>\mu\text{W}</math></td> <td>1.7 <math>\mu\text{W}</math></td> <td>170 nW</td> </tr> <tr> <td><u>Gain setting (high speed) (V/A)</u></td> <td><math>10^3</math></td> <td><math>10^4</math></td> <td><math>10^5</math></td> <td><math>10^6</math></td> <td><math>10^7</math></td> <td><math>10^8</math></td> </tr> <tr> <td>Upper cut-off frequency (-3 dB)</td> <td>175 MHz</td> <td>80 MHz</td> <td>14 MHz</td> <td>3.5 MHz</td> <td>1.8 MHz</td> <td>220 kHz</td> </tr> <tr> <td>Rise/fall time (10 % - 90 %)</td> <td>2.3 ns</td> <td>3.45 ns</td> <td>26 ns</td> <td>94 ns</td> <td>233 ns</td> <td>1.6 <math>\mu\text{s}</math></td> </tr> <tr> <td>NEP (<math>\sqrt{\text{V/Hz}}</math>, @ 850 nm)</td> <td>231 pW</td> <td>10 pW</td> <td>2.2 pW</td> <td>670 fW</td> <td>228 fW</td> <td>76 fW</td> </tr> <tr> <td>Measured at</td> <td>18 MHz</td> <td>8 MHz</td> <td>1.4 MHz</td> <td>350 kHz</td> <td>180 kHz</td> <td>22 kHz</td> </tr> <tr> <td>Integr. input noise (RMS)*</td> <td>4.5 <math>\mu\text{W}</math></td> <td>440 nW</td> <td>31 nW</td> <td>4.8 nW</td> <td>1.3 nW</td> <td>100 pW</td> </tr> <tr> <td>CW saturation power (@ 850 nm)</td> <td>1.7 mW</td> <td>170 <math>\mu\text{W}</math></td> <td>17 <math>\mu\text{W}</math></td> <td>1.7 <math>\mu\text{W}</math></td> <td>170 nW</td> <td>17 nW</td> </tr> </table> <p>* The integrated input noise is measured with a shaded input in the full bandwidth ("FBW") setting (referred to 850 nm). The measurement bandwidth is <math>3 \times</math> the upper cut-off frequency at the specific gain setting; filter slope is a 1st order roll-off.</p> <p>The input referred peak-peak noise can be calculated from the RMS noise as follows:</p> $P_{\text{Input noise peak-to-peak}} = P_{\text{Input noise RMS}} \times 6$ <p>The output noise is given by:</p> $U_{\text{Output noise RMS}} = P_{\text{Input noise RMS}} \times \text{gain}$ $U_{\text{Output noise peak-to-peak}} = U_{\text{Output noise RMS}} \times 6 = P_{\text{Input noise RMS}} \times \text{gain} \times 6$ <p>The integrated noise will be reduced considerably by setting the low pass filter to "1 MHz" or "10 MHz" instead of "FBW". This is especially useful for continuous wave (CW) measurements.</p>		Test conditions	$V_S = \pm 15\text{ V}$ , $T_A = 25\text{ }^\circ\text{C}$ , output load impedance $50\ \Omega$ , warm-up 20 minutes (min. 10 minutes recommended)						Gain	Transimpedance gain	$1 \times 10^2 \dots 1 \times 10^8\ \text{V/A}$ (output load $50\ \Omega$ )					Gain accuracy	$\pm 1\%$ electrical, between settings					Frequency Response	Lower cut-off frequency	DC / 100 Hz, switchable					Upper cut-off frequency (-3 dB)	up to 200 MHz (see table below), switchable to 1 MHz or 10 MHz					Input	Optical CW saturation power	see table below					Noise equivalent power (NEP)	see table below					Performance depending on Gain Setting	<u>Gain setting (low noise) (V/A)</u>	$10^2$	$10^3$	$10^4$	$10^5$	$10^6$	$10^7$	Upper cut-off frequency (-3 dB)	200 MHz	80 MHz	14 MHz	3.5 MHz	1.8 MHz	220 kHz	Rise/fall time (10 % - 90 %)	1.9 ns	3.25 ns	26 ns	92 ns	235 ns	1.6 $\mu\text{s}$	NEP ( $\sqrt{\text{V/Hz}}$ , @850 nm)	322 pW	25 pW	2.9 pW	740 fW	260 fW	78 fW	Measured at	20 MHz	8 MHz	1.4 MHz	350 kHz	180 kHz	22 kHz	Integr. input noise (RMS)*	7.5 $\mu\text{W}$	580 nW	35 nW	4.9 nW	1.3 nW	100 pW	CW saturation power (@ 850 nm)	10 mW	1.7 mW	170 $\mu\text{W}$	17 $\mu\text{W}$	1.7 $\mu\text{W}$	170 nW	<u>Gain setting (high speed) (V/A)</u>	$10^3$	$10^4$	$10^5$	$10^6$	$10^7$	$10^8$	Upper cut-off frequency (-3 dB)	175 MHz	80 MHz	14 MHz	3.5 MHz	1.8 MHz	220 kHz	Rise/fall time (10 % - 90 %)	2.3 ns	3.45 ns	26 ns	94 ns	233 ns	1.6 $\mu\text{s}$	NEP ( $\sqrt{\text{V/Hz}}$ , @ 850 nm)	231 pW	10 pW	2.2 pW	670 fW	228 fW	76 fW	Measured at	18 MHz	8 MHz	1.4 MHz	350 kHz	180 kHz	22 kHz	Integr. input noise (RMS)*	4.5 $\mu\text{W}$	440 nW	31 nW	4.8 nW	1.3 nW	100 pW	CW saturation power (@ 850 nm)	1.7 mW	170 $\mu\text{W}$	17 $\mu\text{W}$	1.7 $\mu\text{W}$	170 nW	17 nW
Test conditions	$V_S = \pm 15\text{ V}$ , $T_A = 25\text{ }^\circ\text{C}$ , output load impedance $50\ \Omega$ , warm-up 20 minutes (min. 10 minutes recommended)																																																																																																																																																		
Gain	Transimpedance gain	$1 \times 10^2 \dots 1 \times 10^8\ \text{V/A}$ (output load $50\ \Omega$ )																																																																																																																																																	
	Gain accuracy	$\pm 1\%$ electrical, between settings																																																																																																																																																	
Frequency Response	Lower cut-off frequency	DC / 100 Hz, switchable																																																																																																																																																	
	Upper cut-off frequency (-3 dB)	up to 200 MHz (see table below), switchable to 1 MHz or 10 MHz																																																																																																																																																	
Input	Optical CW saturation power	see table below																																																																																																																																																	
	Noise equivalent power (NEP)	see table below																																																																																																																																																	
Performance depending on Gain Setting	<u>Gain setting (low noise) (V/A)</u>	$10^2$	$10^3$	$10^4$	$10^5$	$10^6$	$10^7$																																																																																																																																												
	Upper cut-off frequency (-3 dB)	200 MHz	80 MHz	14 MHz	3.5 MHz	1.8 MHz	220 kHz																																																																																																																																												
	Rise/fall time (10 % - 90 %)	1.9 ns	3.25 ns	26 ns	92 ns	235 ns	1.6 $\mu\text{s}$																																																																																																																																												
	NEP ( $\sqrt{\text{V/Hz}}$ , @850 nm)	322 pW	25 pW	2.9 pW	740 fW	260 fW	78 fW																																																																																																																																												
	Measured at	20 MHz	8 MHz	1.4 MHz	350 kHz	180 kHz	22 kHz																																																																																																																																												
	Integr. input noise (RMS)*	7.5 $\mu\text{W}$	580 nW	35 nW	4.9 nW	1.3 nW	100 pW																																																																																																																																												
	CW saturation power (@ 850 nm)	10 mW	1.7 mW	170 $\mu\text{W}$	17 $\mu\text{W}$	1.7 $\mu\text{W}$	170 nW																																																																																																																																												
	<u>Gain setting (high speed) (V/A)</u>	$10^3$	$10^4$	$10^5$	$10^6$	$10^7$	$10^8$																																																																																																																																												
	Upper cut-off frequency (-3 dB)	175 MHz	80 MHz	14 MHz	3.5 MHz	1.8 MHz	220 kHz																																																																																																																																												
	Rise/fall time (10 % - 90 %)	2.3 ns	3.45 ns	26 ns	94 ns	233 ns	1.6 $\mu\text{s}$																																																																																																																																												
	NEP ( $\sqrt{\text{V/Hz}}$ , @ 850 nm)	231 pW	10 pW	2.2 pW	670 fW	228 fW	76 fW																																																																																																																																												
	Measured at	18 MHz	8 MHz	1.4 MHz	350 kHz	180 kHz	22 kHz																																																																																																																																												
	Integr. input noise (RMS)*	4.5 $\mu\text{W}$	440 nW	31 nW	4.8 nW	1.3 nW	100 pW																																																																																																																																												
	CW saturation power (@ 850 nm)	1.7 mW	170 $\mu\text{W}$	17 $\mu\text{W}$	1.7 $\mu\text{W}$	170 nW	17 nW																																																																																																																																												
	Detector	<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 30%;">Detector type</td> <td>Si-PIN photodiode</td> </tr> <tr> <td>Active area</td> <td>1 mm <math>\times</math> 1 mm</td> </tr> <tr> <td>Spectral range</td> <td>400 - 1000 nm</td> </tr> <tr> <td>Sensitivity</td> <td>0.58 A/W typ. (@ 850 nm)</td> </tr> <tr> <td>Dark current</td> <td>0.12 nA typ.</td> </tr> </table>		Detector type	Si-PIN photodiode	Active area	1 mm $\times$ 1 mm	Spectral range	400 - 1000 nm	Sensitivity	0.58 A/W typ. (@ 850 nm)	Dark current	0.12 nA typ.																																																																																																																																						
	Detector type	Si-PIN photodiode																																																																																																																																																	
	Active area	1 mm $\times$ 1 mm																																																																																																																																																	
	Spectral range	400 - 1000 nm																																																																																																																																																	
	Sensitivity	0.58 A/W typ. (@ 850 nm)																																																																																																																																																	
	Dark current	0.12 nA typ.																																																																																																																																																	
Output	<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 30%;">Output voltage rang</td> <td><math>\pm 1\ \text{V}</math> (@ <math>50\ \Omega</math> output load), for linear amplification</td> </tr> <tr> <td>Output impedance</td> <td><math>50\ \Omega</math> (designed for <math>50\ \Omega</math> load)</td> </tr> <tr> <td>Max. output current</td> <td><math>\pm 40\ \text{mA}</math> (short-circuit proof)</td> </tr> <tr> <td>Slew rate</td> <td>1000 V/<math>\mu\text{s}</math></td> </tr> <tr> <td>Output offset compensation</td> <td>adjustable by offset potentiometer and external control voltage, output offset compensation range min. <math>\pm 100\ \text{mV}</math></td> </tr> </table>		Output voltage rang	$\pm 1\ \text{V}$ (@ $50\ \Omega$ output load), for linear amplification	Output impedance	$50\ \Omega$ (designed for $50\ \Omega$ load)	Max. output current	$\pm 40\ \text{mA}$ (short-circuit proof)	Slew rate	1000 V/ $\mu\text{s}$	Output offset compensation	adjustable by offset potentiometer and external control voltage, output offset compensation range min. $\pm 100\ \text{mV}$																																																																																																																																							
Output voltage rang	$\pm 1\ \text{V}$ (@ $50\ \Omega$ output load), for linear amplification																																																																																																																																																		
Output impedance	$50\ \Omega$ (designed for $50\ \Omega$ load)																																																																																																																																																		
Max. output current	$\pm 40\ \text{mA}$ (short-circuit proof)																																																																																																																																																		
Slew rate	1000 V/ $\mu\text{s}$																																																																																																																																																		
Output offset compensation	adjustable by offset potentiometer and external control voltage, output offset compensation range min. $\pm 100\ \text{mV}$																																																																																																																																																		

## 200 MHz Variable Gain Photoreceiver


Specifications (continued)								
DC Monitor Output	Monitor output gain	<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="border-bottom: 1px solid black;">Mode</td> <td style="border-bottom: 1px solid black;">Monitor gain</td> </tr> <tr> <td>Low noise</td> <td>Gain setting divided by -1</td> </tr> <tr> <td>High speed</td> <td>Gain setting divided by -10</td> </tr> </table>	Mode	Monitor gain	Low noise	Gain setting divided by -1	High speed	Gain setting divided by -10
Mode	Monitor gain							
Low noise	Gain setting divided by -1							
High speed	Gain setting divided by -10							
	Monitor output polarity Monitor output voltage range Monitor output bandwidth Monitor output impedance	inverting ±1 V (@ ≥1 MΩ load) DC ... 1 kHz 1 kΩ (designed for ≥1 MΩ load)						
Indicator LED	Function	overload						
Digital Control	Control input voltage range Control input current Overload output	LOW bit: -0.8 V ... +1.2 V, HIGH bit: +2.3 V ... +12 V 0 mA @ 0 V, 1.5 mA @ +5 V, 4.5 mA @ +12 V non active: <0.4 V @ 0 ... -1 mA active: typ. 5 ... 5.1 V @ 0 ... 2 mA						
Ext. Offset Control	Control voltage range Offset control input impedance	±10 V 15 kΩ						
Optical Input Connector	Material FST flange Material FST coupler ring	1.4305 stainless steel, nickel-plated 1.4305 stainless steel, glass bead blasted						
Power Supply	Supply voltage Supply current	±15 V (±14.75 V ... ±16.5 V) ±110 / -90 mA typ. (depends on operating conditions, recommended power supply capability min. ±200 mA)						
Case	Weight Material	360 g (0.79 lbs) AlMg4.5Mn, nickel-plated						
Temperature Range	Storage temperature Operating temperature	-40 °C ... +80 °C 0 °C ... +60 °C						
Absolute Maximum Ratings	Optical input power (CW) Digital control input voltage Analog control input voltage Power supply voltage	12 mW -5 V/+16 V relative to digital ground DGND (pin 9) ±15 V relative to analog ground AGND (pin 3) ±20 V						
Connectors	Input Output Power supply	1.035"-40 threaded flange for free space applications BNC jack (female) LEMO® series 1S, 3-pin fixed socket (mating plug type: FFA.1S.303.CLAC52)						
		 <div style="display: flex; justify-content: space-between; margin-top: 5px;"> <div style="text-align: left;"> <p>PIN 2 -Vs</p> <p>PIN 1 +Vs</p> <p>PIN 3 GND</p> </div> <div style="text-align: right;"> <p>Pin 1: +15 V Pin 2: -15 V Pin 3: GND</p> </div> </div>						

## 200 MHz Variable Gain Photoreceiver

Connectors (continued)

Control port

Sub-D 25-pin, female, qual. class 2



- Pin 1: +12 V (stabilized power supply output\*)
- Pin 2: -12 V (stabilized power supply output\*)
- Pin 3: AGND (analog ground for pins 1 - 8)
- Pin 4: +5 V (stabilized power supply output\*)
- Pin 5: digital output: overload (referred to pin 3)
- Pin 6: DC Monitor output
- Pin 7: NC
- Pin 8: offset control voltage input
- Pin 9: DGND (ground for digital control pins 10 - 16)
- Pin 10: digital control input: gain, LSB
- Pin 11: digital control input: gain
- Pin 12: digital control input: gain, MSB
- Pin 13: digital control input: AC/DC
- Pin 14: digital control input: high speed / low noise
- Pin 15: upper cut-off frequency limit 10 MHz
- Pin 16: upper cut-off frequency limit 1 MHz
- Pin 17 - 25: NC

\*stabilized power supply output current  
 ±12 V: max. ±20 mA, +5V: max. 30 mA

Remote Control Operation

General

Remote control input bits are opto-isolated and connected by a logical OR function to the local switch settings. For remote control set the corresponding local switches to "Remote", "DC", "L" (low noise mode) and "FBW", and select the desired setting via a bit code at the corresponding digital inputs.  
 Mixed operation, e.g. local AC/DC setting and remote controlled gain setting, is also possible.

Gain setting

	Low noise	High speed		
Gain (V/A)	Gain (V/A)	Pin 12	Pin 11	Pin 10
Pin 14=LOW	Pin 14=HIGH	MSB		LSB
10 <sup>2</sup>	10 <sup>3</sup>	LOW	LOW	LOW
10 <sup>3</sup>	10 <sup>4</sup>	LOW	LOW	HIGH
10 <sup>4</sup>	10 <sup>5</sup>	LOW	HIGH	LOW
10 <sup>5</sup>	10 <sup>6</sup>	LOW	HIGH	HIGH
10 <sup>6</sup>	10 <sup>7</sup>	HIGH	LOW	LOW
10 <sup>7</sup>	10 <sup>8</sup>	HIGH	LOW	HIGH

AC/DC setting

Coupling	Pin 13
DC	LOW
AC	HIGH

Low pass filter setting


Upper cut-off freq. limit	Pin 15	Pin 16
full bandwidth	LOW	LOW
10 MHz	HIGH	LOW
1 MHz	LOW	HIGH

High speed / low noise setting

Mode	Pin 14
low noise mode	LOW
high speed mode	HIGH

## 200 MHz Variable Gain Photoreceiver

Scope of Delivery	OE-300-SI-10, internally threaded coupler ring, LEMO® 3-pin connector, datasheet, transport package	
Ordering Information	OE-300-SI-10-FST	1.035"-40 threaded flange for free space applications and for use with various types of optical standard accessories




## 200 MHz Variable Gain Photoreceiver

Typical Performance Characteristic


Frequency response

$$V_{\text{Supply}} = \pm 15 \text{ V}_{\text{DC}}; R_{\text{Load}} = 50 \Omega$$


Gain setting:  $L10^2, H10^3$




Gain setting:  $L10^3, H10^4$




Gain setting:  $L10^4, H10^5$




Gain setting:  $L10^5, H10^6$



Gain setting:  $L10^6, H10^7$



Gain setting:  $L10^7, H10^8$




## 200 MHz Variable Gain Photoreceiver


Typical Performance  
Characteristic (continued)

Input noise equivalent power (NEP)


Gain setting  $L10^2, H10^3$




Gain setting  $L10^3, H10^4$




Gain setting:  $L10^4, H10^5$




Gain setting:  $L10^5, H10^6$



Gain setting:  $L10^6, H10^7$



Gain setting:  $L10^7, H10^8$






## 200 MHz Variable Gain Photoreceiver

Typical Performance  
Characteristic (continued)


### Signal pulse response

Gain setting L10<sup>2</sup>




Rise: 1.85 ns Fall: 1.89 ns

Gain setting H10<sup>3</sup>




Rise: 2.23 ns Fall: 2.27 ns

Gain setting L10<sup>3</sup>




Rise: 3.20 ns Fall: 3.23 ns

Gain setting H10<sup>4</sup>




Rise: 3.44 ns Fall: 3.47 ns

Gain setting L10<sup>4</sup>



Rise: 26.87 ns Fall: 25.66 ns

Gain setting H10<sup>5</sup>




Rise: 27.02 ns Fall: 26.10 ns

## 200 MHz Variable Gain Photoreceiver


Typical Performance  
Characteristic (continued)

Gain setting L10<sup>5</sup>




Rise: 91.80 ns Fall: 91.88 ns

Gain setting H10<sup>6</sup>




Rise: 94.44 ns Fall: 93.16 ns

Gain setting L10<sup>6</sup>




Rise: 233.36 ns Fall: 238.40 ns

Gain setting H10<sup>7</sup>




Rise: 231.92 ns Fall: 234.40 ns

Gain setting L10<sup>7</sup>



Rise: 1606.4 ns Fall: 1584.8 ns

Gain setting H10<sup>8</sup>




Rise: 1621.6 ns Fall: 1608.8 ns

200 MHz Variable Gain Photoreceiver

Dimensions

OE-300-SI-10-FST



all dimensions in mm unless otherwise noted

DZ-OE-300-FST\_R1

FEMTO Messtechnik GmbH  
 Klosterstr. 64  
 10179 Berlin · Germany  
 Phone: +49 30 280 4711-0  
 Fax: +49 30 280 4711-11  
 Email: info@femto.de  
 www.femto.de

Specifications are subject to change without notice. Information provided herein is believed to be accurate and reliable. However, no responsibility is assumed by FEMTO Messtechnik GmbH for its use, nor for any infringement of patents or other rights of third parties which may result from its use. No license is granted by implication or otherwise under any patent or patent rights of FEMTO Messtechnik GmbH. Product names mentioned may also be trademarks used here for identification purposes only.

© by FEMTO Messtechnik GmbH · Printed in Germany